


Dual Flow Contrast Injection for
Coronary CTA Improves
Visualization of the Right Heart

MEDRAD CT

Corey T. Jensen MD
James G. Ravenel MD
Shaun A. Nguyen MD, MPH
Philip Costello MD
U. Joseph Schoepf MD
MUSC, Charleston, USA.




Dual Flow Contrast Injection for Coronary CTA Improves Visualization of the Right Heart

Objective

A saline chaser is ordinarily used at coronary CTA for contrast savings, homogenous attenuation, and reduced streak artifacts. In many cases, however, the void of contrast precludes right heart analysis. We aimed at improving right heart visualization by automated injection of a contrast/saline mixture during the second phase of injection.

Methods

Coronary CTA imaging was performed on 24 patients. Eight(8) patients were scanned with a monophasic, iodine-only protocol using a single-syringe injector; eight(8) patients with a biphasic protocol (iodine bolus followed by a saline chaser) using a dual-syringe injector; and eight(8) patients with a triphasic protocol with DualFlow (MEDRAD) technology with simultaneous injection from two syringes to achieve a desired mixing ratio of contrast and saline. The iodine bolus was followed by a 50:50 saline/contrast ratio and a saline chaser.




Evaluation

Two radiologists rated the visualization of right and left heart structures (heart valves, moderator band, etc.) and the degree of artifacts. One observer performed attenuation measurements of the left and right heart and of the coronary arteries.

Results

Right heart structures were rated significantly better and artifacts occurred less frequently in the DualFlow group, while left heart structures showed no difference. Contrast attenuation in the right heart was significantly lower in the biphasic group than in the monophasic and DualFlow groups. For the coronary arteries, there were no significant differences between the three groups.

Right Heart Structures	Right Heart Contrast Attenuation	Coronary Arteries
DualFlow significantly better with fewer artifacts [P<0.05]	(mean±SE) [P<0.05, One-Way ANOVA] Biphasic group (217.0±69.0) Monophasic group (342.3±37.1) Stellant DualFlow group (322.3±56.2)	No significant difference [P>0.05, Duncan's multiple comparison]



Conclusion

DualFlow injection provides sufficient enhancement for assessment of the right heart while generally avoiding streak artifacts from dense contrast material. Thus, for coronary CTA for reliable detection of right heart pathology (thrombo-emboli, tumors, etc.), the triphasic injection approach with Stellant D's DualFlow appears recommendable.

MEDRAD reserves the right to modify the specifications and features described herein, or discontinue manufacture of the product described at any time without prior notice or obligation. Please contact your authorized MEDRAD representative for the most current information.

MEDRAD, and Stellant are federally registered trademarks, and Performance. For life. is a trademark of MEDRAD, INC. U.S.A.



The Malcolm Baldrige National Quality Award is our nation's premier award for performance excellence and quality achievement. Established in 1987 and named after the former Secretary of Commerce, the Baldrige Award recognizes exemplary achievements in seven areas: leadership, strategic planning, customer and market focus, measurement, analysis and knowledge management, human resource focus, process management and business results. The U.S. Commerce Department's National Institute of Standards and Technology (NIST) manages the Baldrige National Quality Program in close cooperation with the private sector. The award is traditionally presented by the President of the United States in a special Washington, D.C. ceremony.

MEDRAD, INC.
One Medrad Drive
Indianola, PA 15051-0780 USA
412-767-2400

Customer Service/Orders
1-800-MEDRAD-1
(1-800-633-7231)

Customer Service FAX
(412) 767-4120

International FAX
(412) 767-4128

For Service Repairs
1-800-MEDRAD-7
(1-800-633-7237)

MEDRAD Service FAX
(412) 767-4126

www.medrad.com/stellant
Email address for more
information:
info@medrad.com



© 2006 MEDRAD, INC. All Rights Reserved.

CJ CS EN 204109 Rev A



MEDRAD®
Performance. For life.™