



As of January 2016

	MEDRAD® Mark 7 Arterion (with VFlow)	MEDRAD® Avanta
Installation Options	<ul style="list-style-type: none"> > Pedestal > Floor Mount > Table Mount > Overhead Counterpoise 	<ul style="list-style-type: none"> > Pedestal > Table Mount
Injectable Fluid	<ul style="list-style-type: none"> > X-ray contrast media 	<ul style="list-style-type: none"> > X-ray contrast media > Saline
Syringe Size	<ul style="list-style-type: none"> > 150 mL > Single patient use only 	<ul style="list-style-type: none"> > 150 mL > Up to 5 consecutive patient use
Max. Pressure Limit	<ul style="list-style-type: none"> > 1,200 psi 	<ul style="list-style-type: none"> > 1,200 psi
Contrast Media Flow Rate (Fixed)	<ul style="list-style-type: none"> > 0.1–45.0 mL/s in 0.1 ml/sec increments > 0.1–59.9 mL/min in 0.1 ml/min increments 	<ul style="list-style-type: none"> > 0.1–45.0 mL/s in 0.1 ml/sec increments
Contrast Media Flow Rate (Variable)	<ul style="list-style-type: none"> > 1–10 mL/s in 0.1 ml/sec increments > V-Flow option required 	<ul style="list-style-type: none"> > 1–10 mL/s in 0.1 ml/sec incremental
Saline Flush Capability	<ul style="list-style-type: none"> > – 	<ul style="list-style-type: none"> > Available
Saline Flow Rate (Fixed)	<ul style="list-style-type: none"> > – 	<ul style="list-style-type: none"> > 0.75–1.50 mL/s
Injection Volume	<ul style="list-style-type: none"> > 1~150 mL 	<ul style="list-style-type: none"> > 1~150 mL
Dual Display Option	<ul style="list-style-type: none"> > Available 	<ul style="list-style-type: none"> > –
Protocol Memory	<ul style="list-style-type: none"> > 40 	<ul style="list-style-type: none"> > 40
Injection History Record	<ul style="list-style-type: none"> > 50 	<ul style="list-style-type: none"> > N/A > Can be obtained by Bayer Engineer

MEDRAD® Mark 7 Arterion
Injection System

MEDRAD® Avanta
Fluid Management Injection System



Clear Direction. From Diagnosis to Care.

Every day it's your expertise that provides clear direction on a patient's clinical journey. Putting an end to uncertainty. We recognize the difference you make. That's why we're committed to supporting you in your crucial task. Every single day.

More information on radiology.bayer.com

All rights reserved. This publication or parts thereof may not be translated into other languages or reproduced in any form mechanical or electronic (including photocopying, tape recording, microcopying) or stored in a data carrier or computer system without written consent of Bayer Pharma AG. Bayer reserves the right to modify the specifications and features described herein, or discontinue manufacture of the product described at any time without prior notice or obligation. Please contact your authorized Bayer representative for the most current information. Bayer, the Bayer Cross, and MEDRAD are trademarks of the Bayer group of companies. © 2016 Bayer



Bayer

Bayer Pharma AG
13342 Berlin
Germany

More information on
radiology.bayer.com

UNITED STATES
Bayer Medical Care Inc.
1 Bayer Drive
Indianola, PA 15051-0780
U.S.A.
Phone: +1.412.767.2400
+1.800.633.7231
Fax: +1.412.767.4120

EUROPE
Bayer Medical Care B.V.
Horsterweg 24
6199 AC Maastricht Airport
The Netherlands
Phone: +31 (0) 43-3583600
Fax: +31 (0) 43-3656598